

LiteStar®

Disposable anesthesia face masks.



D-288-2010

LiteStar anesthesia face mask is for disposable use. The teardrop shape of the mask supports an anatomical fit. The transparent mask body allows for continuous observation of the patient's condition. There is a color coded reinforcing ring placed around the mask connector, which provides stable connection with breathing circuit systems (sizes 3-6).

The "Traffic Light" color coded reinforcing rings are for easier size identification. The smaller size (3) has red color, mid-sizes (4-5) have yellowish color and the large size (6) has green color.

ADDITIONAL MASK ACCESSORIES

Additional hook rings and headgears for fixation are available.

Easy application is supported through additional mask accessories: hook rings and headgears for fixation. Three different headgears in small and large designs offer patients and physicians a soft and elastic fixation.




D-957-2010

LiteStar® – Disposable anesthesia face masks.

ORDER INFORMATION

Description	Package of	Order No.
Mask, LiteStar, without hook ring		
size 1, neonate	30 pcs.	MP01501
size 2, infant	30 pcs.	MP01502
size 3, toddler	30 pcs.	MP01503
size 4, small adult	30 pcs.	MP01504
size 5, medium adult	30 pcs.	MP01505
size 6, large adult	30 pcs.	MP01506
Headgear, Silicone, reusable		
small	1 pc.	MP01557
large	1 pc.	MP01558
Headgear, Soft Fabric, disposable		
small	1 pc.	MP01559
large	1 pc.	MP01560
Headgear, Soft Rubber, disposable		
small	1 pc.	MP01561
large	1 pc.	MP01562
Hook ring for LiteStar masks		
with 15 mm OD connector (size 1–2)	30 pcs.	MP01551
with 22 mm ID connector (size 3–6)	30 pcs.	MP01552



HEADQUARTERS
Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

**REGION EUROPE CENTRAL
AND EUROPE NORTH**
Dräger Medical GmbH
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH
Dräger Médical S.A.S.
Parc de Haute
Technologie d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
dlmfr-contact@draeger.com

**REGION MIDDLE EAST, AFRICA,
CENTRAL AND SOUTH AMERICA**
Dräger Medical GmbH
Branch Office Dubai
Dubai Healthcare City, P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 436 24 762
Fax +971 436 24 761
contactuae@draeger.com

REGION ASIA / PACIFIC
Dräger Medical
South East Asia Pte Ltd
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION NORTH AMERICA
Dräger Medical, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 215 721 5400
Toll-free +1 800 437 2437
Fax +1 215 723 5935
info.usa@draeger.com

Manufacturer:
Dräger Medical GmbH
23542 Lübeck, Germany
The quality management system at
Dräger Medical GmbH is certified
according to ISO 13485, ISO 9001
and Annex II.3 of Directive 93/42/EEC
(Medical devices).