

Quickstop hinged arms

Keeping everything close at hand. To leave your hands free to care for the patient, all the joints of the Quickstop hinged arm can be fixed into place and released with just one quick hand-movement.



The hinged arm is totally flexible, allowing it to be positioned precisely where it is needed. All cables, lines and hose systems can be arranged in such a way as to more or less completely prevent accidental disconnection, making the Quickstop hinged arm an essential accessory in anaesthesia, intensive care, surgery or indeed just about any other area of the hospital. The arm is provided with a universal rail bracket which can be used to attach it to all kinds of different objects (operating tables, beds, wall-mounted rails

and devices). Thanks to the high quality of its materials, you will appreciate the benefits of this durable and maintenance-free Dräger product for years to come. Quickstop hinged arms - keeping everything close at hand.

- Quick
- Precise
- Flexible
- Universal
- Robust



Quickstop hinged arms

TECHNICAL DATA


Quickstop hinged arm 1	Order no.
Length 1067 mm, contains: 3 x hose clamps (2 pcs.) 1 x hose clamps, coaxial (2 pcs.) Load capacity: max. 2 kg	2M 85 705
Quickstop hinged arm 2	Order no.
Length 1078 mm, contains: 2 x hose clamps (2 pcs.) 1 x hose clamps, coaxial (double) 1 x hose clamp with ball Load capacity: max. 2 kg Optional accessories Holder for Y-piece (for ball mounting)	2M 85 706 84 09 746
Quickstop hinged arm Pod	Order no.
Length 460 mm For the fixation of Dräger monitoring pods. Load capacity: max. 1,5 kg Optional accessories Quickstop set for the fixation of MultiMed and NeoMed pods.	2M 86 464 MP 00 653
Quickstop hinged arm, short	Order no.
Length 403 mm, contains: 1 x universal tube holder plate for several cable- and tube dimensions Load capacity: max. 4 kg	2M 86 463
Spare Parts	
Hose clamps, coaxial (double)	MP00493
Hose clamp (double)	8406138
Hose clamp with ball	8409841
Holder for Y-piece	8409746
Universal tube holder plate (for ball mounting)	2M86465
Quickstop set for the fixation of MultiMed and NeoMed pods.	MP00653



Quickstop hinged arm 1



Quickstop hinged arm 2



Holder for Y-piece




Quickstop hinged arm Pod




Quickstop hinged arm, short



Quickstop hinged arm, short



MT-0814-2008



MT-0814-2008

HEADQUARTERS
Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23558 Lübeck, Germany

www.draeger.com

REGION EUROPE CENTRAL AND EUROPE NORTH
Dräger Medical GmbH
Moislinger Allee 53-55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH
Dräger Médical S.A.S.
Parc de Haute
Technologie d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
dlmfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA, CENTRAL AND SOUTH AMERICA
Dräger Medical GmbH
Branch Office Dubai
Dubai Healthcare City, P.O. Box 505108
Dubai, United Arab Emirates
Tel + 971 436 24 762
Fax + 971 436 24 761
contactuae@draeger.com

REGION ASIA / PACIFIC
Dräger Medical
South East Asia Pte Ltd
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION NORTH AMERICA
Draeger Medical, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 215 721 5400
Toll-free +1 800 437 2437
Fax +1 215 723 5935
info.usa@draeger.com

Manufacturer:
Dräger Medical GmbH
23542 Lübeck, Germany
The quality management system at Dräger Medical GmbH is certified according to ISO 13485, ISO 9001 and Annex II.3 of Directive 93/42/EEC (Medical devices).