

## Caddy and CompactCaddy Transportation solution for Oxylog® 1000 and 2000

Saving a person's life in an emergency depends on your expertise and the devices you work with. Practical transportation solutions are essential – for you and your equipment.



### FREE YOUR HANDS

The Caddy from Dräger Medical offers you an ideal solution for transporting the Oxylog 1000 and Oxylog 2000. The Caddy's simple and smooth design means you can carry, transport and use the Oxylog comfortably and efficiently. At the same time, the Caddy protects your equipment in rough conditions. So much convenience integrated into one simple frame – less juggling of equipment, more time to care for your patients!

### SELECTING THE TRANSPORTATION SOLUTION YOU PREFER

Since you have to carry the load, you are the best person to judge what suits you. That's why we offer two versions of the Caddy:

### Caddy

For a complete, stand-alone ventilation system that will hold the Oxylog, all its accessories and an oxygen cylinder, choose the Caddy. As the cylinder protection bar is adjustable in length and reversible, you can be sure your cylinders will fit – even if you prefer to mount in the opposite direction.

### CompactCaddy

If you prefer a multi-functional use for your cylinders and want to carry the oxygen cylinder separately, choose the Compact-Caddy.



An ideal solution for transporting the Oxylog® 1000 und Oxylog® 2000.



#### CHOOSING VERSATILE VENTILATION

The Caddy enables you to choose and use the best Oxylog ventilator for your mission. Mount just one Dräger wall holder in your ambulance and you can quickly and securely attach the Oxylog models 1000 and 2000. Why compromise when one solution meets all your needs?

#### ENJOYING HANDS-FREE TRANSPORTATION

A belt is provided so you can carry the Caddy on your shoulder or back. But you don't always have to carry the Caddy – its claw is designed so it can be hung on a variety of rails, tubes, beds and stretchers.



MT-3414-2003

MT-9233-2009



**PROTECTING YOU  
AND YOUR EQUIPMENT**

The Caddy frame is a robust construction designed to withstand tough conditions. The surfaces are smooth and rounded for your safety; the optimized point of gravity and fixed handle ensure you can comfortably carry an Oxylog, cylinder and all the accessories you need. When not in use, the Caddy's protective bag can be easily pulled over the entire system.

**Profiting from a cost-effective solution.**

Both Caddy models offer a complete, versatile and cost-effective transportation solution:

- Ergonomic design to optimize carrying your load – an important health and safety consideration
- Versatile construction suitable for multiple cylinder sizes – no change of frame needed
- Quick and secure mounting onto the modular wall holder – saves time, reduces risk and fits Oxylog 1000 and Oxylog 2000 models
- Easy attachment to various rails, poles, beds and stretchers
- Longer-lasting equipment thanks to sturdy steel frame and protection bag



A belt is provided so you can carry the Caddy on your shoulder or back.

MT-3382-2003

**TECHNICAL DATA****W x H x D**

Caddy, incl. Oxylog 1000	510 x 430 x 250 mm / 20.1 x 16.9 x 9.84 inch
Caddy, incl. Oxylog 2000	510 x 430 x 280 mm / 20.1 x 16.9 x 11.0 inch
CompactCaddy, incl. Oxylog 1000	280 x 430 x 280 mm / 11.0 x 16.9 x 11.0 inch
CompactCaddy, incl. Oxylog 2000	280 x 350 x 230 mm / 11.0 x 13.8 x 9.06 inch

**Weight**

Caddy	Approx. 3.8 Kg / 8.38 lbs
CompactCaddy	Approx. 2.0 Kg / 4.41 lbs

**Suitable for cylinder sizes (fits D-size cylinder):**

Diameter	80 -120 mm/ 3.15 - 4.72 inch
Length (incl. valve and pressure reducer)	400 to 580 mm/ 15.7 - 22.8 inc

Shock and vibration	In accordance with EN 794-3, EN 1789, EN 13718 (vehicles)
Classification (MDD93/42/EEC)	Class I
UMDNS Code1	11-835
All-round Wall Holder (W x H x D)	300 x 460 x 85 mm / 11.8 x 18.1 x 3.35 inch

**NOTE:**

The Oxylog 2000 itself was discontinued in 2008, but the Caddy and Compact Caddy can still be purchased as a spare part.

**HEADQUARTERS**

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53-55  
23558 Lübeck, Germany

[www.draeger.com](http://www.draeger.com)

**REGION EUROPE CENTRAL  
AND EUROPE NORTH**

Dräger Medical GmbH  
Moislinger Allee 53-55  
23558 Lübeck, Germany  
Tel +49 451 882 0  
Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

**REGION EUROPE SOUTH**

Dräger Médical S.A.S.  
Parc de Haute  
Technologie d'Antony 2  
25, rue Georges Besse  
92182 Antony Cedex, France  
Tel +33 1 46 11 56 00  
Fax +33 1 40 96 97 20  
[dlmfr-contact@draeger.com](mailto:dlmfr-contact@draeger.com)

**REGION MIDDLE EAST, AFRICA,  
CENTRAL AND SOUTH AMERICA**

Dräger Medical GmbH  
Branch Office Dubai  
Dubai Healthcare City, P.O. Box 505108  
Dubai, United Arab Emirates  
Tel + 971 436 24 762  
Fax + 971 436 24 761  
[contactuae@draeger.com](mailto:contactuae@draeger.com)

**REGION ASIA / PACIFIC**

Dräger Medical  
South East Asia Pte Ltd  
25 International Business Park  
#04-27/29 German Centre  
Singapore 609916, Singapore  
Tel +65 6572 4388  
Fax +65 6572 4399  
[asia.pacific@draeger.com](mailto:asia.pacific@draeger.com)

**REGION NORTH AMERICA**

Dräger Medical, Inc.  
3135 Quarry Road  
Telford, PA 18969-1042, USA  
Tel +1 215 721 5400  
Toll-free +1 800 437 2437  
Fax +1 215 723 5935  
[info.usa@draeger.com](mailto:info.usa@draeger.com)

**Manufacturer:**

Dräger Medical GmbH  
23542 Lübeck, Germany  
The quality management system at  
Dräger Medical GmbH is certified  
according to ISO 13485, ISO 9001  
and Annex II.3 of Directive 93/42/EEC  
(Medical devices).